



Oxygen Tank

Dose

Your doctor has prescribed you a flow rate of _____liters.

Using an Oxygen Tank

1. Place the regulator on a full tank.
2. Check for leaks around the regulator when changing tanks.
3. Check the pressure gauge to ensure adequate oxygen supply. Change the tanks when the pressure is below 500 psi.
4. Dial in the liter flow. You are prescribed_____ liters.
5. Keep tanks lying on the floor or secured in an appropriate cart or case.
6. Check to make sure there is flow running through the nasal cannula.

Continuous Flow Regulator

Your regulator is a continuous flow regulator which provides you with a steady flow of oxygen. This flow can be adjusted by dialing it on your regulator.

