



Oxlife Liberty 2 Portable Oxygen Concentrator

Dose

Your doctor has prescribed you a flow rate of _____.

What is a Portable Oxygen Concentrator?

The Liberty 2 Portable Oxygen Concentrator is small, compact, lightweight, and portable. It works by taking the air we breathe, filtering out the nitrogen, and making it concentrated oxygen. The oxygen is then delivered by a pulse dose or continuously through a nasal cannula by the machine.

What is a pulse dose?

A pulse dose delivers a burst of oxygen every time you inhale. This will allow you to get oxygen at your own pace. The Liberty 2 can deliver a flow of 1–10-liter pulse doses.

What is a continuous dose?

A continuous dose of oxygen gives you a constant flow of oxygen all the time. The Liberty 2 is designed to give a continuous flow of oxygen ranging from 0.5 – 2 liters.

Continuous Flow Settings

| Flow Setting | Liter Flow |
|---------------------|-------------------|
| Setting 1 | 0.5 LPM |
| Setting 2 | 0.75 LPM |
| Setting 3 | 1 LPM |
| Setting 4 | 1.25 LPM |
| Setting 5 | 1.5 LPM |
| Setting 6 | 1.75 LPM |
| Setting 7 | 2 LPM |

How to switch between continuous and pulse doses and change flow of oxygen

To switch between continuous and pulse mode simply hold down the “M” button. The mode will appear on the screen. Use the “+” to increase the flow and “-“ button to decrease the flow of the Liberty 2.

Battery Life

Battery life is dependent on settings of oxygen. See charts below for details.

Pulse Flow Battery Life

| Flow Setting | Hours of Use |
|--------------|------------------|
| 1 | Up to 4.5 Hours |
| 2 | Up to 4 Hours |
| 3 | Up to 3.5 Hours |
| 4 | Up to 3 Hours |
| 5 | Up to 2.5 Hours |
| 6 | Up to 2 Hours |
| 7 | Up to 1.75 Hours |
| 8 | Up to 1.75 Hours |
| 9 | Up to 1.25 Hours |
| 10 | Up to 1.25 Hours |

Continuous Flow Battery Life

| Flow Setting | Hours of Use |
|----------------------|------------------|
| Setting 1 (0.5 LPM) | Up to 3.25 Hours |
| Setting 2 (0.75 LPM) | Up to 2.5 Hours |
| Setting 3 (1 LPM) | Up to 2 Hours |
| Setting 4 (1.25 LPM) | Up to 1.75 Hours |
| Setting 5 (1.5 LPM) | Up to 1.5 Hours |
| Setting 6 (1.75 LPM) | Up to 1.25 Hours |
| Setting 7 (2 LPM) | Up to 1 Hour |

You will be provided with a car charger and a wall charger. When you are away from home, please keep your concentrator on charge as much as you can and be mindful of the battery life displayed on the screen. It takes 2.5 hours to fully charge a battery. You can use the Liberty 2 when it is charging.

Cleaning

Clean the outside of the case with a mild soap and water. Do NOT let any liquids into the controls, interior of the case, and oxygen tubing.

Filters allow for adequate air flow. Clean the external filter weekly to ensure this. Use a mild soap and water. Make sure the external filter is completely dry before putting it back on your concentrator.



